SFWRENG 3K04: Software Development

Assignment 2 – Part 2 – DCM Design

Dr. Alan Wassyng

December 4, 2019

Group 5: More Life Pacemaker

Adam Bujak – bujaka – 400113347

Eric Hillebrand – hillebre – 400143468

Harneet Singh – singhh76 – 400110275

Karan Gill – gillk20 – 400130973

 $Stephen\ Scott - scotts 24 - 400139933$

Table of Contents

Likely Requirement Changes	
Likely Design Changes	3
Modules	4
Application	4
GUI	4
GUI Controller	4
GUI Abstraction Layer	7
Data Management	10
User Account Manager	10
Database Manager	14
SQLite ORM Library	17
SQLite Database	17
Pacemaker Communication	17
Protocol	17
DCM Communication Controller	17
DCM Serial Manager	20
Appendix	23
Module Flowchart	23

Likely Requirement Changes

- The user interface shall be capable of visually indicating when a different pacemaker device is approached than was previously interrogated
- The Set Clock function shall set the data and time of the device
- The New Patient function shall allow a new device to be interrogated without exiting the software application
- Upon the user's request:
 - o A Bradycardia Parameters Report shall be available
 - A Temporary Parameters Report shall be available

Likely Design Changes

- The current GUI has some issues where certain features aren't very intuitive. For instance, if all fields are blank the user can still press program, and the last saved data in the database will be programmed to the pacemaker. This can be confusing and may be rectified in later versions.
- Add additional comments onto the interface in order to make it more user friendly
- Change the DCM interface such that it is easier to distinguish between the programmable parameters and the echoed parameters
- Implement a way for the DCM to detect which COM port the pacemaker is on

Modules

Application

Description

This is the main application module for all other components of the Device Controller-Monitor (DCM). It starts the Graphical User Interface (GUI) and handles the transmission of data between modules.

GUI

GUI Controller

Description

The GUI controller handles the drawing of the custom interface to the display by calling the GUI Library through the GUI Abstraction Layer (see Appendix for a flowchart of the modules). The GUI controller takes care of tasks such as starting menus, filling forms with values, and building the GUI.

Public Functions

Class	Functions/Methods
LoginMenuData()	init(self, userNameLabel, passwordLabel, signInButtonText, newUserButtonText)
	setCallbacks(self, callbacks)
ProgramMenuData()	init(self, fieldLabels, buttonTexts, dropDownLabelText, dropDownOptions, currentOption)
	setCallbacks(self, callbacks)
CreateUserData()	init(self, userNameLabel, passwordLabel, createUserButtonText, cancelButtonText)
	setCallbacks(self, callbacks)
Screen	init(self, screenname, data)
GUIC	init(self) setCallbacks(self, callbacks) startGUI(self) drawScreen(self, screen) updateGUI(self) getLoginData(self) getNewUserData(self) getProgramData(self) drawFirstScreen(self) p_drawLoginScreen(self, data) p_drawProgrammingScreen(self, data) p_drawCreateUserScreen(self, data)

Black Box Behaviour

Function	Input(s)	Output(s)
ALLinit() FUNCTIONS	Parameters for a class structure	Sets the values of the structure with the given inputs. Returns nothing
ALL setCallbacks(self, callbacks) FUNCTIONS	Callbacks	Enables button functionality/interaction. Returns nothing
startGUI(self)	None	Draws the first screen. Returns nothing
drawScreen(self, screen)	Screen class	Draws a screen of either type Login, Programming, or Create User. Returns nothing
updateGUI(self)	None	Updates the interface. Returns nothing
getLoginData(self)	Input field data	Returns GUI field data on Login Screen in a LoginData class
getNewUserData(self)	Input field data	Returns GUI field data on Create User Screen in a LoginData class
getProgramData(self)	Input field data	Returns GUI field data on Programming Screen in a ProgrammedData class
p_drawFirstScreen(self)	None	Draws the Login screen. Returns nothing
p_drawLoginScreen(self, data)	Data/text to fill the Login screen	Draws the Login screen. Returns nothing
<pre>p_drawProgrammingScreen(self, data)</pre>	Data/text to fill Programming screen	Draws the Programming screen. Returns nothing
<pre>p_drawCreateUserScreen(self, data)</pre>	Data/text to fill Create User screen	Draws the Create User screen. Returns nothing

Global Variables

Variable(s)	Description		
LoginData	Imported class from the User Account Manager module, each of which contains a username and password for a user		
ApplicationCallbacks	Imported class from a file that contains numerous button callbacks		
Enum	Imported Enum from enum		
loginScreen	A Screen class created with ScreenNames(0) and the LoginMenuData class		
programmingScreen	A Screen class created with ScreenNames(1) and the ProgramMenuData class.		
createUserMenuScreen	A Screen class created with ScreenNames(2) and the CreateUserData class.		

Function	Description
ALLinit() FUNCTIONS	For all classes this is a constructor function that is automatically called with the creation of every instance of a class. This function allows the class to initialize itself with the attributes/values it is given.
ALL setCallbacks(self, callbacks) FUNCTIONS	Tkinter's button widget provides a command callback when a user clicks a button. This function enables the functionality for all relevant classes without the need to create extra GUI objects.
startGUI(self)	Uses the p_drawFirstScreen method to start the GUI
drawScreen(self, screen)	Clears the window, determines which screen type should be drawn (LOGIN_SCREEN, PROGRAMMING_SCREEN, or CREATE_USER_SCREEN), and then uses the appropriate method to draw the screen.
updateGUI(self)	Updates the GUI
getLoginData(self)	Retrieves data from the GUI input fields when the Login Screen is up and returns the data in a LoginData class
getNewUserData(self)	Retrieves data from the GUI input fields when the Create User Screen is up and returns the data in a LoginData class
getProgramData(self)	Retrieves data from the GUI input fields when the Programming Screen is up and returns the data in a ProgrammedData class
p_drawFirstScreen(self)	Draws the first screen (Login Screen) of the GUI by using self.drawScreen(loginScreen)
p_drawLoginScreen(self, data)	Draws the Login Screen with field labels, button texts, and button callbacks.
<pre>p_drawProgrammingScreen(self, data)</pre>	Draws the Programming Screen with drop-down labels, drop-down options, field labels, button texts, and button callbacks.
<pre>p_drawCreateUserScreen(self, data)</pre>	Draws the Create User Screen with field labels, button texts, and button callbacks.

GUI Abstraction Layer

Description

The GUI Abstraction Layer hides the specific working details of drawing the GUI from the GUI controller in order to create a Separation of Concerns; this helps allow for the reuse and independent maintenance/upgrades of modules.

Public Functions

Class	Functions/Methods
GUIAL	init(self)
	update(self)
	setTitle(self, title)
	clearWindow(self)
	drawTwoFieldsTwoButtonLayout(self, fieldLabels, buttonTexts, buttonCallbacks)
	drawNFieldsNButtonsOneDropDownLayout(self, dropDownLabelText, currentDropDownItem, dropDownOptions, fieldLabels, buttonTexts, buttonCallbacks)
	getEntryData(self)
	displayLabelEntry(self, programMode, dropDownOptions, fieldLabels) getProgramMode(self)
	displayErrorMessageLoginS(self, errorCode)
	displayErrorMessageProgramS(self, errorCodeRate, errorCodeChamber)
	setNEntryData(self, data)

Function	Input(s)	Output(s)
init()	Parameters for a class	Sets the values of the
	structure	structure with the
		given inputs. Returns
		nothing
update(self)	None	Updates the GUI.
		Returns nothing
setTitle(self, title)	Title string	Writes the input in the
		title bar of the GUI

Function	Input(s)	Output(s)
clearWindow(self)	None	Clears the window. Returns nothing
drawTwoFieldsTwoButtonLayout(self, fieldLabels, buttonTexts, buttonCallbacks)	Input field text, button text, and button callbacks	Draws a screen with two user input fields and buttons, all with associated text. Returns nothing
drawNFieldsNButtonsOneDropDownLayout(self, dropDownLabelText, currentDropDownItem, dropDownOptions, fieldLabels, buttonTexts, buttonCallbacks)	Drop down menu label, value of the current drop-down option, the set of drop-down options, an array of input field texts, button text, and button callbacks	Draws a screen with ten user input fields, a drop-down menu, and buttons, all with associated texts. Returns nothing
displayLabelEntry(self, programMode, dropDownOptions, fieldLabels)	The program's current state, the options for the drop-down menu, and the text labels to display on the input fields	Sets the label entries to their respective values. Returns nothing
getProgramMode(self)	None	Returns the program's current state
displayErrorMessageLoginS(self, errorCode)	The relevant error code for the login screen	Displays a login error message
displayErrorMessageProgramS(self, errorCodeRate, errorCodeChamber)	The relevant error code parameters	Displays a programming error message
getEntryData(self)	None	Returns an array of the data inputted into the data fields of the current screen
setNEntryData(self, data)	Array of values to be set	Sets each value to the entries in the data array. Returns nothing

Variable(s)		Description
Variables	From	The file dcm_constants.py contains the variables that hold the text to
src.dcm_constants		be displayed when drawing various screens (i.e.
		C_LOGIN_USERNAME_LABEL = "Username",
		C_PROGRAM_BUTTON_TEXT = "Program")

Function	Description
init()	Constructor function that is automatically called with the creation of every instance of a class. This function allows the class to initialize itself with the attributes/values it is given
update(self)	Updates the GUI and initializes it if it hasn't been already
setTitle(self, title)	Sets the text to be shown in the title bar of the GUI if it has been initialized
clearWindow(self)	Clears all elements on the current screen by parsing through the list that is the output of the getEntryData(self) function
drawTwoFieldsTwoButtonLayout(self, fieldLabels, buttonTexts, buttonCallbacks)	Draws two user input fields and a button to the screen. fieldLabels is an array of labels to be used for the input fields. buttonTexts is an array of texts to be displayed in the button. buttonCallbacks is an array of button callback functions. These are drawn using the GUI Library (tkinter)
drawNFieldsNButtonsOneDropDownLayout(self, dropDownLabelText, currentDropDownItem, dropDownOptions, fieldLabels, buttonTexts, buttonCallbacks)	Draws ten user input fields and a button to the screen. dropDownLabelText is the label for dropDownMenu. currentDropDownItem is the value of the current dropDownOption. dropDownOptions is the set of dropDownOptions. fieldLabels is an array of labels to be used for the input fields. buttonTexts is an array of texts to be displayed in the button. buttonCallbacks is an array of button callback functions. These are drawn using the GUI Library (tkinter). This method creates a tkinter variable and uses it to set the currentDropDownItem. It also uses a tkinter link function to change the dropdown.
displayLabelEntry(self, programMode, dropDownOptions, fieldLabels)	Displays the necessary label entry by comparing the inputs with the program mode. The functions goes through a series of "if" statements where the programMode is compared with the values in the dropDownOptions[] array.
getProgramMode(self)	Returns the program's current mode
displayErrorMessageLoginS(self, errorCode)	Sets the label to be displayed with the error code text inside it. Returns nothing
displayErrorMessageProgramS(self,	Sets the labels to be displayed with the error code
errorCodeRate, errorCodeChamber)	texts within. Returns nothing
getEntryData(self)	Gets the information inputted into the data field and returns it in an array
setNEntryData(self, data)	Sets the values in the parameter array to be single strings that the user inputs using tk.Entry

Data Management

User Account Manager

Description

The User Account Manager handles the creation of and logging in/out of Device Control-Monitor users. It also sets the devices programmable parameters: the upper/lower rate limits and the atrium/ventricle pulse amplitudes, pulse widths, sensing thresholds, and refractory periods.

Public Functions

Class	Functions/Methods	
LoginData	init(self, p_username, p_password)	
ProgammedData	init(self, p_upperRateLim, p_lowerRateLim, p_atriumPulseAmp, p_atriumPulseWidth, p_atriumSensThres, p_atriumRefracPeriod, p_ventriclePulseAmp, p_ventriclePulseWidth, p_ventricleRefracPeriod) p_lowerRateLim, p_lowerRateLim, p_atriumSensThres, p_ventriclePulseAmp, p_ventriclePulseAmp, p_ventricleSensThres, p_ventricleRefracPeriod)	
DUAM	init(self) getSessionState(self) signInUser(self, p_loginData) signOut(self) p_makeAdminUser(self) makeNewUser(self, p_loginData, p_adminPassword) changeUserPassword(self, p_username, p_existingPassword, p_newPassword) validUser(self) validNumUsers(self) programRateLim(self, p_upperRateLim, p_lowerRateLim) programAtriaPara(self, p_atriumPulseAmp, p_atriumPulseWidth, p_atriumSensThres, p_atriumRefracPeriod) programVentriclePara(self, p_ventriclePulseAmp, p_ventriclePulseWidth, p_ventriclePulseWidth, p_ventriclePares, p_ventricleSensThres, p_ventricleRefracPeriod)	
None	hash_password(password) verify_password(stored_password, provided_password)	
None	frange(start, stop, step)	

Function	Input(s)	Output(s)
ALLinit() FUNCTIONS	Parameters for a class structure	Sets the values of the structure with the given inputs. Returns nothing
getSessionState(self)	None	Returns the state of the session (logged in/out)
signInUser(self, p_loginData)	Username and password from input fields	If the credentials are valid, the user is logged in and a success code is returned. Else, a failure code of "invalid credentials" is returned if either the username or password are invalid
signOut(self)	None	Signs the user out of the DCM. Returns True
p_makeAdminUser(self)	None	Creates the admin user in the database. Returns nothing
makeNewUser(self, p_loginData, p_adminPassword)	Login data (username & password) and admin password	If successful, it creates a user in the database with a username, password, and role and a success code is returned. Else, a failure code is returned denoting "missing permissions", "too many users", or "existing user"
changeUserPassword(self, p_username. P_existingPassword, p_newPassword)	Username, current password, and the new password	If successful, changes the old password of a user with the new one given as input. Else, returns a failure code denoting "invalid credentials" if either the username or password are incorrect
validUser(self)	None	If logged out, returns False. Else, returns True
validNumUsers(self)	None	Returns False if the number of users in the database is greater than or equal to 10. Else, returns True
programRateLim(self, p_upperRateLim, p_lowerRateLim)	Input field data for the upper and lower rate limits	Sets the user's data with the appropriate inputs. Returns nothing

Function	Input(s)	Output(s)
programAtriaPara(self,	Input field data for the atrium's	Sets the user's data with
p_atriumPulseAmp,	pulse amplitude, pulse width,	the appropriate inputs.
p_atriumPulseWidth,	sensing threshold, and refractory	Returns nothing
p_atriumSensThres,	period	
p_atriumRefracPeriod)		
programVentriclePara(self,	Input field data for the ventricle's	Sets the user's data with
p_ventriclePulseAmp,	pulse amplitude, pulse width,	the appropriate inputs.
p_ventriclePulseWidth,	sensing threshold, and refractory	Returns nothing
p_ventricleSensThres,	period	
p_ventricleRefracPeriod)		
hash_password(password)	Password	Returns the hashed
		password
verify_password(stored_password,	A password in the database and a	Returns True if the inputs
provided_password)	password from the user	are the same. Else, False
frange(start, stop, step)	3 float numbers	Returns a float number

Variable(s)		Description
Variables	From	The imported file dcm_constants.py contains the variables that hold
src.dcm_constants		the text to be displayed when drawing various screens (i.e. C_LOGIN_USERNAME_LABEL = "Username", C_PROGRAM_BUTTON_TEXT = "Program")
Variables DCMDatabase.dbpm	From	This imported file contains the User class which hold the username, password, and userRole. It also has the DBPM class (Database Peewee Manager)

Function	Description
ALLinit() FUNCTIONS	For all classes this is a constructor function that is automatically called with the creation of every instance of a class. This function allows the class to initialize itself with the attributes/values it is given. The init function for the DUAM class however also creates the admin user when the database is empty,
	and the admin user has never been created yet.
getSessionState(self)	Returns the current state of the user (logged out/in)
signInUser(self, p_loginData)	Signs the user in. It checks that the user exists in the database and then fetches the user data to compare with the inputted data and decide whether or not to log the user in.
signOut(self)	Signs the user out
p_makeAdminUser(self)	Adds the admin user to the database if it is not already there

Function	Description
makeNewUser(self, p_loginData, p_adminPassword)	Adds a new user and its login credentials to the database. Returns a failure code if there are too many users, if the user already exists, or if the user is invalid
changeUserPassword(self, p_username. P_existingPassword, p_newPassword)	Changes a user's stored password through the database manager after checking that the user exists and verifying the existing password
validUser(self)	Checks if the current user is valid by returning False if signed out or if the user object is None
validNumUsers(self)	Checks if the number of users in the data base is greater than or equal to 10. Returns False if this is the case
programRateLim(self, p_upperRateLim, p_lowerRateLim)	Sets the current user's upper and lower rate limits in the database given the inputs
programAtriaPara(self, p_atriumPulseAmp, p_atriumPulseWidth, p_atriumRefracPeriod) p_atriumSensThres,	Sets the current user's atrium pulse amplitude, pulse width, sensing threshold, and refractory period in the database given the inputs
programVentriclePara(self, p_ventriclePulseAmp, p_ventriclePulseWidth, p_ventricleSensThres, p_ventricleRefracPeriod)	Sets the current user's ventricle pulse amplitude, pulse width, sensing threshold, and refractory period in the database given the inputs
hash_password(password)	Hashes a password for storing
verify_password(stored_password,	Verifies a stored password against one provided by
provided_password)	the user
frange(start, stop, step)	Performs the same function as the regular Range() function but with floating point numbers. If the inputs "stop" and "step" are null, then "start" is set to 0.0 and "step" to 1.0

Database Manager

Description

The Database Manager handles the actual working details of storing, retrieving, and changing user information by communicating with the database.

Public Functions

Class	Functions/Methods
User	init(self, p_username, p_password, p_userRole) getUsername(self) getPassword(self) getRole(self)
UserProgramData	init(self, p_upperRateLim, p_atriumPulseAmp, p_atriumPulseAmp, p_atriumPulseWidth, p_atriumRefracPeriod, p_ventriclePulseWidth, p_ventriclePulseWidth, p_ventricleRefracPeriod) There are 10 functions for getting all of these values and another 10 for setting these values after initialization. For simplicity's sake, they will be listed next to each other. getLowerRateLimit(self) / setLowerRateLimit(self) getAtrialAmplitude(self) / setAtrialAmplitude(self) getAtrialPulseWidth(self) / setAtrialPulseWidth(self) getAtrialSensingThreshold(self) / setAtrialRefractoryPeriod(self) getVentricularAmplitude(self) / setVentricularAmplitude(self) getVentricularPulseWidth(self) / setVentricularIPulseWidth(self) getVentricularSensingThreshold(self) / setVentricularIPulseWidth(self) getVentricularSensingThreshold(self) / setVentricularIPulseWidth(self) getVentricularSensingThreshold(self) / setVentricularIPulseWidth(self) getVentricularSensingThreshold(self) / setVentricularRefractoryPeriod(self) getVentricularRefractoryPeriod(self) / setVentricularRefractoryPeriod(self)

Class	Functions/Methods
DBPM	init(self)
	closeDatabase(self)
	getDatabaseInstance(self)
	createUser(self, p_username, p_password, p_role)
	userExists(self, p_username)
	getUserData(self, p_username)
	getNumUsers(self)

Private Functions

	Class	Functions/Methods
DBPM		p_createDataTables(self)

Function	Input(s)	Outputs(s)
ALLinit() FUNCTIONS	Parameters for the class structure	Sets the values of the structure with the given inputs. Returns nothing
getUsername(self)	None	Returns the username of the user
getPassword(self)	None	Returns the password of the user
getRole(self)	None	Returns the role of the user (doctor, nurse, admin, etc.)
closeDatabase(self)	None	Returns with the method that closes the database
getDatabaseInstance(self)	None	Returns with the method that gets the database instance
createUser(self, p_username, p_password, p_role, p_data)	Username, password, and a role	Creates the user in the database and stores the data in the database. Returns nothing
userExists(self, p_username)	Username	Returns True if a username is the database matches the input. Else, returns False
getUserData(self, p_username)	Username	Returns a User class with username, password, and role if a username matching the input is found. Else, returns None
getNumUsers(self)	None	Returns the number of users in the database

Function	Input(s)	Outputs(s)
p_createDataTables(self)	None	Creates the tables used in the application. Returns nothing
For the set functions mentioned in the Public Functions table	None	Returns the value corresponding to the function name
For the get functions mentioned in the Public Functions table	User input from a data field	Sets the value from the input field corresponding to the field label and the function name. Returns nothing

Function		Input(s)
Variables	From	The imported file dcm_constants.py contains the variables that hold
src.dcm_constants		the text to be displayed when drawing various screens (i.e.
		C_LOGIN_USERNAME_LABEL = "Username",
		C_PROGRAM_BUTTON_TEXT = "Program")
database		SqliteDatabase structure with a given path, size, and keys

Function	Description
ALLinit(self) FUNCTIONS	For all classes this is a constructor function that is automatically called with the creation of every instance of a class. This function allows the class to initialize itself with the attributes/values it is given. The init function for the DBMP class also connects to the database and creates data tables.
getUsername(self)	Returns the username of the user
getPassword(self)	Returns the password of the user
getRole(self)	Returns the role of the user (doctor, nurse, admin, etc.)
closeDatabase(self)	Closes the database
getDatabaseInstance(self)	Returns the current database instance
createUser(self, p_username, p_password, p_role, p_data)	Creates a new user and data variable in the database. It stores the data from the input variable into their appropriate locations as well as the username, password, and role of the new user
userExists(self, p_username)	Checks if a user exists by looping through users in the database and comparing usernames. Returns True or False
getUserData(self, p_username)	Loops through user data in the database looking for a username matching the input. If one is found, returns a User class (with username, password, role, and data)

Function	Description
getNumUsers(self)	Loops through user data in the database and counts the number of users
p_createDataTables(self)	Creates all tables used in the application

SQLite ORM Library

An ORM is an Object-Relational-Mapper. It acts as a layer between the SQLite database and the object-oriented code (database manager). The one used for this DCM is called Peewee and works with both Python and SQLite. This allows us to use the database without having to transform objects into a suitable format for it.

SQLite Database

SQLite is a software library that provides a relational database management system (RDBMS). RDBMSes are a type of database manager that store data in a format using rows and columns.

Pacemaker Communication

Protocol

Description

The standard communication protocol is to send a start byte, followed by a command byte, followed by data of the form '3B2H2B5H2B' where B represents an unsigned uint8 and H represents a uint16.

Example

To program the pacemaker parameters the DCM would send a start byte (0x16), then the pacemaker parameter programming command byte (0x55), then 21 bytes of parameter data.

DCM Communication Controller

Description

This module acts as an abstraction layer for the DCM Serial Manager. It handles the communication between the DCM and the Pacemaker by opening/closing then connection and reading/writing the values stored on the Pacemaker.

Public Functions

Class	Functions/Methods

DCC	init(self)
	deinit(self)
	openConnection(self)
	closeConnection(self)
	programPacemaker(self, params)
	getPacemakerData(self)
	getElectrogram(self)

Private Functions

Class	Eurostiana/Mathada
Class	Functions/Methods

0.20.55	
DCC	p_prependDataWithStartCode(self, data)
	p_transmitData(self, data)
	p_convertToInts(self, number, numberOfBytes)
	p_convertPacemakerParamsToByteArray(self, pacemakerParams)
	p_sendEchoCommand(self)

Function	Input(s)	Outputs(s)
init(self)	None	Initializes serialManager. Returns nothing
deinit(self)	None	Deinitializes serialManager. Returns nothing
openConnection(self)	None	Initializes serialManager. Returns nothing
closeConnection(self)	None	Deinitializes serialManager. Returns nothing
p_prependDataWithStartCode(self, data)	The data to be sent to the Pacemaker	Returns the data prepended with the serial start byte
p_transmitData(self, data)	The data to be sent to the Pacemaker	Sends data to the Pacemaker and prints a fail/success statement. Statements correspond with the possible returns of True, False, or nothing.

Function	Input(s)	Outputs(s)
programPacemaker(self, params)	Programmable parameters for the Pacemaker	Converts params into the byte array. Can return False
p_convertToInts(self, number, numberOfBytes)	Number to be converted and the number of bytes	Returns the converted number in an array.
p_convertPacemakerParamsToByteArray(self, pacemakerParams)	Programmable parameters for the Pacemaker	Returns a byte array of the Pacemaker parameters
p_sendEchoCommand(self)	None	Returns the success of the data transmit; True or False
getPacemakerData(self)	None	Returns the parameters store in the Pacemaker in a byte array
getElectrogram(self)	None	Returns an array of electrogram values from the Pacemaker

Variable(s)		Description
Variables	From	The imported file dcm_constants.py contains the variables that hold
src.dcm_constants		the text to be displayed when drawing various screens (i.e.
		C_LOGIN_USERNAME_LABEL = "Username",
		C_PROGRAM_BUTTON_TEXT = "Program")
Variables	from	The imported file datatypes.py contains general variables used
common.datatypes		throughout the DCM. Most of these are the Pacemaker parameters
		and the various programming modes the Pacemaker can be in.

Function	Description		
init(self)	Initializes serialManager.		
deinit(self)	Deinitializes serialManager.		
openConnection(self)	Opens the serial connection by initializing the serialManager		
closeConnection(self)	Closes the serial connection by deinitializing the serialManager		
p_prependDataWithStartCode(self, data)	Prepends the serial communication start byte to the byte array. Returns that array		
p_transmitData(self, data)	Prepares the Pacemaker data using p_prependDataWithStartCode(self, data) then sends it using the DCM Serial Manager		
programPacemaker(self, params)	Sends the Pacemaker programmable parameters after converting them into a byte array		
p_convertToInts(self, number, numberOfBytes)	Converts the inputted number using bit shifting and logical operations		

Function	Description	
p_convertPacemakerParamsToByteArray(self, pacemakerParams)	Converts the Pacemaker parameters into a byte array by using p_convertToInts(self, number, numberOfBytes) and adding the returns from each function call to the array	
p_sendEchoCommand(self)	Sends an echo command to the Pacemaker using p_transmitData(self, data)	
getPacemakerData(self)	Gets the parameters stored on the Pacemaker in a byte array using the serialManager	
getElectrogram(self)	Reads the electrogram from the Pacemaker and returns the values in an array	

DCM Serial Manager

Description

This module handles the lower level details and the minute organization of communication necessary for the Pacemaker to send and receive data from the DCM.

Public Functions

Class	Functions/Methods
DSM	init(self)
	deinit(self)
	write(self, data)
	writeString(self, dataStr)
	readUntil(self, expected)
	readLine(self)
	getSerialPort(self)
	checkSerialPort(self)

Function	Input(s)	Outputs(s)
init(self)	None	Initializes serialManager. Can print "Error Opening Com Port" if an error occurs. Returns nothing
deinit(self)	None	Deinitializes the serialManager. Returns nothing

Function	Input(s)	Outputs(s)
write(self, data)	Input byte array to write	Prints the serial data and its length. Writes this data to the serial port. Can return a failure code
writeString(self, dataStr)	Input string to write	Writes the byte array to the serial port. Returns nothing
readUntil(self, expected)	Expected end character	Reads the serial data until the expected character. Returns this data; can also return a failure code
read (self, n)	Number of bytes to read	Reads n number of bytes of serial data. Returns this data; can also return a failure code.
readLine(self)	None	Reads and returns one line of serial data. Can return a failure code
getSerialPort(self)	None	Returns the name of the serial port. Can return a failure code
checkSerialPort(self)	None	Returns true if the serial port is available. Else, returns a failure code

Variable(s)		Description	
Variables src.dcm_constants	From	The imported file dcm_constants.py contains the variables that hold the text to be displayed when drawing various screens (i.e.	
		C_LOGIN_USERNAME_LABEL = "Username", C_PROGRAM_BUTTON_TEXT = "Program")	
Variables common.failCodes	from	The imported file failCodes.py contains the enumeration for the possible errors that may occur such as there being too many users, invalid rate input, etc	

Function	Description
init(self)	Initializes the serialManager or prints an error message if it cannot for some reason
deinit(self)	Deinitializes the serialManager by closing the serial port
write(self, data)	Writes the data to the serial port and prints it along with its length to the screen. Returns a failure code if the port is not open/available
writeString(self, dataStr)	Writes the data string to the serial port using write(self, data)

Function	Description
readUntil(self, expected)	Reads the serial data until the expected character and returns that data. Returns a failure code if the serial port is not open/available
read (self, n)	Reads n number of bytes of serial data. Returns a failure code if the serial port is not open/available
readLine(self)	Reads one line of serial data using readUntil(self, expected) with no input. Returns a failure code if the serial port is not open/available
getSerialPort(self)	Returns the name of the serial port. Returns a failure code if the serial port is not open/available
checkSerialPort(self)	Checks to see if the serial port is open and available. Returns a failure code if the serial port is not open/available

Appendix

Module Flowchart

